

ST-analyzer

Windows installation

Ver. 0.1.0

Committee

Wonpil Im
Chair, University of Kansas

Jong Cheol Jeong
Developer, University of Kansas

Sunhwan Jo
University of Chicago

Yifei Qi
University of Kansas

All rights reserved @ ImLab

This work is supported by NSF ABI



Objective

This tutorial guides how users can install ST-analyzer into Windows.

Installing ST-analyzer into Windows

Installing ST-analyzer into Windows system requires Linux emulator that makes Linux-like system environment. Therefore, besides installing the emulator, rest of procedures is similar to those in Linux.

1. Install Cygwin,

- a) download setup-x86.exe from <http://www.cygwin.com/>
- b) run setup-x86.exe, select packages: gcc, gfortran, python, python setuptools, lapack and numpy

2. Install scipy

- a) download source code of scipy, blas and atlas, and move them to somewhere in cygwin, e.g, `~/temp`. After copying the files from Windows to Cygwin, you may need to change the permission using `chmod`
- b) unzip blas and atlas
 - `tar -xjvf atlas3.11.11.tar.bz2`
 - `tar -xzvf blas.tgz`
- c) set up environmental variables
 - `export BLAS_SRC=~/temp/BLAS`
 - `export ATLAS_SRC=~/temp/ATLAS`
- d) unzip and install scipy
 - `tar -xzvf scipy-0.12.0.tar.gz`
 - `cd scipy-0.12.0`
 - `python setup.py install`
- **Note:** If you want to compile BLAS and ATLAS by yourself to get optimal performance, check this page:
 - <http://www.scipy.org/scipylib/building/index.html>

3. Install MDAnalysis

- `tar -xzvf MDAnalysis-0.7.7.tar.gz`
- `cd MDAnalysis-0.7.7`
- `python setup.py install`

4. Run ST-Analyzer

- `cd ST-analyzer-master/stanalyzer`
- `python manage.py runserver 8000`